

Invitrogen Superscript II Reverse Transcriptase 1st strand cDNA

- Primer

Pd(N)₆ oligos (e.g. from Invitrogen): stock comes 3 µg/µl
dilute 1/30 in DEPC water to get 100 ng/µl working concentration
20 µl working concentration are enough for a 100 µl RT

The following protocol is written for a 100 µl RT 1st strand cDNA reaction

- In screw-cap tube

20 µl random primers
30 µl RNA (e.g. directly from Amersham QuickPrep Micro mRNA Purification Kit)
4 µl dNTPs (25 mM each)
66 µl DEPC H₂O

□ 65 °C for 5 min, then quickly chill on ice

- collect drops by briefly spinning (6000 rpm quickspin)

add:

40 µl First-Strand Buffer (comes with the kit)
20 µl 0,1 M DTT
10 µl RNAse Inhibitor*

□ mix contents gently, 25°C for 2 min

- add 10 µl (200u/µl) Superscript Transcriptase mix by gently pipetting up&down

□ 25°C for 10 min
□ 42°C for 50 min
□ 70°C for 15 min

* e.g. Invitrogen RNAseOUT or Amersham RNAGuard